

Anti-CCS antibody (1-274) (STJ28886)

STJ28886

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CCS (1-274) is suitable for use in Western Blot.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9973
Gene Symbol	CCS
Uniprot ID	CCS_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-274 of human CCS (NP_005116.1).
Immunogen Region	1-274
Specificity	
Immunogen Sequence	