

Anti-CIDEA antibody (20-219) (STJ29969)

STJ29969

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CIDEA (20-219) is suitable for use in Immunofluorescence.
Applications	IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF 1:50-1:100
Formulation	PBS containing 0.02% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Storage Instruction	Store in a freezer at -20°C and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1149
Gene Symbol	CIDEA
Uniprot ID	CIDEA_HUMAN
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-219 of human CIDEA (NP_001270.1).
Immunogen Region	20-219
Specificity	
Immunogen Sequence	