

Anti-HCG free beta antibody [5012] (STJ400015)

STJ400015

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HCG free beta is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5012
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	37 mM citrate, 125 mM phosphate, pH 6.0, 0.9 % NaCl, 0.095 % NaN3 as a preservative
Isotype	IgG1
Storage	36 months from manufacturing at 2–8°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21 days OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human chorionic gonadotropin free beta subunit
Immunogen Region	
Specificity	
Immunogen Sequence	