

Anti-FABP3 antibody [2304] (STJ400204)

STJ400204

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-FABP3 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2304
Concentration	5.0 mg/ml (+/-10 %)
Conjugation	Unconjugated
Purification	Greater than 95% purity
Dilution Range	
Formulation	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
Isotype	IgG1
Storage	24 months from manufacturing at 2–8°C. Temperature analysis:-70 °C, 21 days N/D-20 °C, 21 days OK +4 °C, 21 days OK +30 °C, 21
Instruction	OK +35 °C, 21 days OK +45 °C, 7 days OK. Shelf life stability testing at 2–8 °C in the product buffer.

TARGET INFORMATION

Gene ID	2170
Gene Symbol	FABP3
Uniprot ID	FABPH_HUMAN
Immunogen	Human fatty acid binding protein 3 (FABP3)
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	