

Anti-FGFR2 antibody (520-570) {Biotin} (STJ501030)

STJ501030

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FGFR2 (520-570) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.75 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10,000 IP: 1:150 IHC: 1:50-1:100
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
Uniprot ID	FGFR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 520-570 on human Fibroblast Growth Factor receptor 2 protein.
Immunogen Region	520-570
Specificity	
Immunogen Sequence	