

Anti-Hepatitis C virus Core NS3 + NS4 antibody (N-Term) {Biotin} (STJ501318)

STJ501318

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hepatitis C virus Core NS3 + NS4 (N-Term) is suitable for use in Dot Blot, ELISA and Western Blot research applications.
Applications	DB, ELISA, WB
Host/Source	Rabbit
Reactivity	Hepatitis C

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.50-0.65 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000-1:50, 000 ELISA: 1:20, 000-1:50, 000
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	N-epitope specific synthetic peptide corresponding to unique amino acid sequences on hepatitis C polyprotein.
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	