

Anti-Outermembrane protein 26 H. pylori antibody {Biotin} (STJ501362)

STJ501362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Outermembrane protein 26 H. pylori is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IP, IHC, WB
Host/Source	Rabbit
Reactivity	H.pylori

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.80-1.35 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:5, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:300 CFM: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Specificity	
Immunogen Sequence	