

Anti-OPTN antibody {Biotin} (STJ501368)

STJ501368

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OPTN is suitable for use in Confocal Microscopy, ELISA, Immunocytochemistry, Immunofluorescence, Immunoprecipitation, Immunohistochemistry and Western Blot research applications.
Applications	CM, ELISA, ICC, IF, IP, IHC, WB
Host/Source	Rabbit
Reactivity	H.pylori

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.80-1.35 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:5, 000 IP: 1:150 IHC: 1:200 ICC: 1:200 IF: 1:300 CFM: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10133
Gene Symbol	OPTN
Uniprot ID	OPTN_HUMAN
Immunogen	Heat killed and denatured Helicobacter pylori extracts
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	