

Anti-MINPP1 antibody (1-100) {Biotin} (STJ501757)

STJ501757

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MINPP1 (1-100) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Human, Mouse, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.51-0.62 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9562
Gene Symbol	MINPP1
Uniprot ID	MINP1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1-100 on rat MINPP1 protein.
Immunogen Region	1-100
Specificity	
Immunogen Sequence	