

Anti-Cell Surface Glycoprotein MUC18 antibody (330-353) {Biotin} (STJ501820)

STJ501820

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cell Surface Glycoprotein MUC18 (330-353) is suitable for use in ELISA, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IP, WB
Host/Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5-1.5 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	78967
Gene Symbol	Mcam
Uniprot ID	MUC18_RAT
Immunogen	Synthetic peptide taken within amino acid region 330-353 on rat Cell Surface Glycoprotein MUC18 protein.
Immunogen Region	330-353
Specificity	
Immunogen Sequence	