

Anti-PDE6B antibody {FITC} (STJ502362)

STJ502362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PDE6B is suitable for use in ELISA, Confocal Microscopy, Immunocytochemistry, Immunofluorescence, Immunohistochemistry and Western Blot research applications.
Applications	ELISA, CM, ICC, IF, IHC, WB
Host/Source	Rabbit
Reactivity	Bovine, Human, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5158
Gene Symbol	PDE6B
Uniprot ID	PDE6B_HUMAN
Immunogen	Synthetic cyclic peptide common to all PDE6 Beta variants
Immunogen Region	
Immunogen Specificity	
Immunogen Sequence	