

Anti-DNA polymerase eta antibody (200-300) {Biotin} (STJ502748)

STJ502748

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNA polymerase eta (200-300) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	S.cerevisiae

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.70-0.98 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 200-300 on Chain A, Yeast DNA polymerase Eta.
Immunogen Region	200-300
Specificity	
Immunogen Sequence	