

## Anti-TICAM2 antibody (140-170) {FITC} (STJ503268)

STJ503268

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-TICAM2 (140-170) is suitable for use in ELISA and Western Blot research applications. |
| <b>Applications</b>      | ELISA, WB   |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Cat, Human, Monkey, Pig   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 0.60-0.65 µg/µl   |
| <b>Conjugation</b>         | FITC  |
| <b>Purification</b>        | Affinity Purified   |
| <b>Dilution Range</b>      | WB: 1:250<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200     |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">100302736/353376</a>   |
| <b>Gene Symbol</b>        | <a href="#">TICAM2</a>   |
| <b>Uniprot ID</b>         | <a href="#">TCAM2_HUMAN</a>  |
| <b>Immunogen</b>          | Synthetic peptide taken within amino acid region 140-170 on human TIRAP3b long form protein. |
| <b>Immunogen Region</b>   | 140-170  |
| <b>Specificity</b>        |  |
| <b>Immunogen Sequence</b> |  |