

Anti-TLR12 antibody (440-470) {FITC} (STJ503293)

STJ503293

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR12 (440-470) is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Rabbit
Reactivity	Chimpanzee, Human, Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.55-0.69 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000
Formulation	
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	384059
Gene Symbol	Tlr12
Uniprot ID	TLR12_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 440-470 from mouse TLR12 protein.
Immunogen Region	440-470
Specificity	
Immunogen Sequence	