

Anti-TLR6 antibody (195-225) {FITC} (STJ503304)

STJ503304

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TLR6 (195-225) is suitable for use in ELISA, Immunohistochemistry, Immunoprecipitation and Western Blot research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	Rabbit
Reactivity	Cat, Human, Mouse, Pig, Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.78-0.80 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:20, 000 ELISA: 1:20, 000 IP: 1:200 IHC: 1:100
Formulation	
Isotype	IgG
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10333
Gene Symbol	TLR6
Uniprot ID	TLR6_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 195-225 on human TLR6 protein.
Immunogen Region	195-225
Specificity	
Immunogen Sequence	