

FEL D1 Major allergen 1 polypeptide chain 1 protein (STJ504065)

STJ504065

GENERAL INFORMATION

Product Type	Proteins
Short Description	FEL D1 Major allergen 1 polypeptide chain 1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 1-50 on feline Major Cat Allergen I polypeptide chain 1.
Immunogen Region	1-50
Specificity	
Immunogen Sequence	