

FEL D1 Major allergen 1 polypeptide chain 2 protein (STJ504066)

STJ504066

GENERAL INFORMATION

Product Type	Proteins
Short Description	FEL D1 Major allergen 1 polypeptide chain 2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken within amino acid region 42-88 on feline Major Cat Allergen I polypeptide chain 1.
Immunogen Region	42-88
Specificity	
Immunogen Sequence	