

SIRT7 protein (STJ504936)
STJ504936**GENERAL INFORMATION**

Product Type	Proteins
Short Description	SIRT7-protein if for use in research applications.
Applications	ELISA, IHC, IP, WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51547
Gene Symbol	SIRT7
Uniprot ID	SIRT7_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 340-390 on human SIRT7 protein.
Immunogen Region	340-390
Specificity	
Immunogen Sequence	