

**TM6SF2 protein (STJ505013)**  
STJ505013**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	TM6SF2-protein if for use in research applications.
<b>Applications</b>	DB, ELISA, IHC, IP, WB
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">107770</a>
<b>Gene Symbol</b>	<a href="#">Tm6sf2</a>
<b>Uniprot ID</b>	<a href="#">TM6S2_MOUSE</a>
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 220-270 on mouse TM6SF2 protein.
<b>Immunogen Region</b>	220-270
<b>Specificity</b>	
<b>Immunogen Sequence</b>	