

C1ORF115 protein (STJ505863)
STJ505863**GENERAL INFORMATION**

Product Type	Proteins
Short Description	C1ORF115-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	Mouse, Rat

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	226777
Gene Symbol	NA
Uniprot ID	CA115_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 50-100 on mouse Chromosome 1 Open Reading Frame 115.
Immunogen Region	50-100
Specificity	
Immunogen Sequence	