

C1ORF115 protein (STJ505876)
STJ505876**GENERAL INFORMATION**

Product Type	Proteins
Short Description	C1ORF115-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	Human

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79762
Gene Symbol	C1orf115
Uniprot ID	CA115_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 75-125 on human Chromosome 1 Open Reading Frame 115.
Immunogen Region	75-125
Specificity	
Immunogen Sequence	