

Anti-Cd19 antibody [1D3] {FITC} (STJA0000004)

STJA0000004

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1D3
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12478
Gene Symbol	Cd19
Uniprot ID	CD19_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse splenocytes were stained with Anti-CD19 antibody (STJA0000004) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.