

Anti-CD11c/Itgax antibody [N418] {APC} (STJA0000064)

STJA0000064

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Integrin alpha-X is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N418
Concentration	
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
Uniprot ID	ITAX_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine splenocytes were stained with Anti-CD11c antibody (STJA0000064) and Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (Left). Splenocytes stained with Anti-Mouse MHC II (I-A/I-E) Monoclonal Antibody (FITC Conjugated) (Right) were used as control

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081