

## Anti-Cd80 antibody [16-10A1] {PE/Cyanine5.5} (STJA0000084)

STJA0000084

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Armenian Hamster monoclonal antibody anti-T-lymphocyte activation antigen CD80 is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC  |
| <b>Host/Source</b>       | Armenian Hamster  |
| <b>Reactivity</b>        | Mouse   |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | 16-10A1  |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | PE/Cyanine5.5  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |            |
|---------------------------|------------|
| <b>Gene ID</b>            | 12519      |
| <b>Gene Symbol</b>        | Cd80       |
| <b>Uniprot ID</b>         | CD80_MOUSE |
| <b>Immunogen</b>          |            |
| <b>Immunogen Region</b>   |            |
| <b>Specificity</b>        |            |
| <b>Immunogen Sequence</b> |            |



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD80 antibody (STJA0000084) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081