

Anti-CD127/II7r antibody [A7R34] {PE/Cyanine5} (STJA0000257)

STJA0000257

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-7 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	A7R34
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16197
Gene Symbol	Il7r
Uniprot ID	IL7RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Alexa Fluor 488 anti-mouse CD3 Alexa Fluor 488 anti-mouse CD3
 C57BL/6 murine splenocytes were stained with Anti-Mouse CD127 Monoclonal Antibody (PE/Cyanine5 Conjugated) (Used at 0.2 μ g/ 10^6 cells dilution) and Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated). Splenocytes stained with Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (right) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081