

Anti-CD11b antibody [M1/70] {PerCP} (STJA0000424)

STJA0000424

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse, Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Itgam
Uniprot ID	ITAM_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b Monoclonal Antibody (PerCP Conjugated) (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.