

Anti-Cd8a antibody [53-6.7] {PerCP} (STJA0000599)

STJA0000599

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-T-cell surface glycoprotein CD8 alpha chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	53-6.7
Concentration	0.2 mg/mL
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12525
Gene Symbol	Cd8a
Uniprot ID	CD8A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Mouse splenocytes were stained with Anti-Mouse CD8a Monoclonal Antibody (PerCP Conjugated) [Used at 0.05 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081