

Anti-Cd4 antibody [OX-38] {PerCP} (STJA0000612)

STJA0000612

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OX-38
Concentration	0.2 mg/mL
Conjugation	PerCP
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24932
Gene Symbol	Cd4
Uniprot ID	CD4_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



Rat splenocytes were stained with Anti-Rat CD4 (domain 1 Conjugated) Monoclonal Antibody (PerCP Conjugated) (Used at 0.1 µg/10⁶ cells dilution) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081