

Anti-CD183/Cxcr3 antibody [CXCR3-173] {AF488} (STJA0000686)

STJA0000686

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Armenian Hamster monoclonal antibody anti-C-X-C chemokine receptor type 3 is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Armenian Hamster |
| Reactivity | Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | CXCR3-173 |
| Concentration | 0.5 mg/mL |
| Conjugation | AF488 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|---------------------------|-------------|
| Gene ID | 12766 |
| Gene Symbol | Cxcr3 |
| Uniprot ID | CXCR3_MOUSE |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



C57BL/6 murine splenocytes were stained with Anti-CD183/CXCR3 antibody (STJA0000686) [Used at 0.2 µg/10⁶ cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081