

Anti-CD41/Itga2b antibody [MWRReg30] {PE/Cyanine7} (STJA0001014)

STJA0001014

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin alpha-IIb is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	MWRReg30
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16399
Gene Symbol	Itga2b
Uniprot ID	ITA2B_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	



C57BL/6 murine platelets were stained with Anti-CD41 antibody (STJA0001014) (filled gray histogram). Unstained platelets (empty black histogram) were used as control.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081