

Anti-CD11c/Itgax antibody [N418] {Biotin} (STJA0001224)

STJA0001224

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Integrin alpha-X is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N418
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
Uniprot ID	ITAX_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	