

Anti-CD86 antibody [BU63] (STJA0001270)

STJA0001270

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-lymphocyte activation antigen CD86 is suitable for use in Flow Cytometry, Western Blot, Universal and Immunoprecipitation research applications.
Applications	FC, WB, Block, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	BU63
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG1k
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	942
Gene Symbol	CD86
Uniprot ID	CD86_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	