

Anti-CD11b antibody [M1/70] (STJA0001523)

STJA0001523

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD11b is suitable for use in Universal, Flow Cytometry, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Neut, FC, IHC, IP
Host/Source	Rat
Reactivity	Human, Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M1/70
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	NA
Gene Symbol	Itgam
Uniprot ID	ITAM_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	