

Anti-CD103/Itgae antibody [M290] (STJA0001550)

STJA0001550

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin alpha-E is suitable for use in Flow Cytometry, Immunohistochemistry and Universal research applications.
Applications	FC, IHC, Block
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M290
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16407
Gene Symbol	Itgae
Uniprot ID	ITAE_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	