

Anti-CD3 Epsilon antibody [145-2C11] (STJA0001660)

STJA0001660

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-T-cell surface glycoprotein CD3 epsilon chain is suitable for use in Flow Cytometry, Western Blot, Immunohistochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	FC, WB, IHC, IF, IP
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	145-2C11
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	12501
Gene Symbol	Cd3e
Uniprot ID	CD3E_MOUSE
Immunogen	
Region	
Specificity	
Immunogen Sequence	