

Anti-HLA-DRA antibody [L243] (STJA0001707)

STJA0001707

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-HLA class II histocompatibility antigen, DR alpha chain is suitable for use in Universal, Universal, Universal, Flow Cytometry, Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Cell Sep-Neg, Block, Binding, FC, WB, IHC, IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	L243
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	3122
Gene Symbol	HLA-DRA
Uniprot ID	DRA_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	