

Anti-CD49b/Itga2 antibody [DX5] {Biotin} (STJA0001741)

STJA0001741

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Integrin alpha-2 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	DX5
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgMk
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID	16398
Gene Symbol	Itga2
Uniprot ID	ITA2_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	