

Anti-CD45.2/Ptprc antibody [104.2] (STJA0001761)

STJA0001761

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Universal, Flow Cytometry and Immunohistochemistry research applications.
Applications	Block, FC, IHC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	104.2
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2ak
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	19264
Gene Symbol	Ptprc
Uniprot ID	PTPRC_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	