

Anti-TCR Beta antibody [H57-597 (HB218)] {Biotin} (STJA0001770)

STJA0001770

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-TCR Beta is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	H57-597 (HB218)
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing at 2-8°C. Do not freeze this product.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence