

Anti-IL2 antibody [JES6-5H4] (STJA0002142)

STJA0002142

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 is suitable for use in Universal, Universal, IConfocal MicroscopyC, Immunohistochemistry and Immunoprecipitation research applications.
Applications	Neut, Stim, ICFC, IHC, IP
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	JES6-5H4
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	
Dilution Range	
Formulation	0.2 µm filtered in PBS, pH7.2. Azide Free (AF) /Low Endotoxin (LE) : Contains no stabilizers or preservatives. Endotoxin level is < 2 EU/mg as determined by LAL gel clotting assay.
Isotype	IgG2bk
Storage	Recommend storing at 2-8°C. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	16183
Gene Symbol	IL2
Uniprot ID	IL2_MOUSE
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	