

Anti-IL-5 antibody [TRFK5] (STJA0002164)

STJA0002164

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rat monoclonal antibody anti-IL-5 is suitable for use in Confocal Microscopy, ELISA, Western Blot and Immunohistochemistry research applications. |
| Applications | ICFC, ELISA, WB, IHC |
| Host/Source | Rat |
| Reactivity | Human, Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Monoclonal |
| Clone ID | TRFK5 |
| Concentration | 0.5 mg/mL |
| Conjugation | Unconjugated |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300 |
| Isotype | IgG1k |
| Storage | Recommend storing at 2-8°C. Do not freeze this product. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|---------------------------|
| Gene ID | 3567 |
| Gene Symbol | IL5 |
| Uniprot ID | IL5_HUMAN |
| Immunogen | |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |