

## Anti-HIV-1 Rev antibody (34-40aa) [10] (STJA0002844)

STJA0002844

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-HIV-1 Rev (34-40aa) is suitable for use in ELISA, Western Blot, Immunofluorescence, Immunoprecipitation and Immunohistochemistry research applications.
<b>Applications</b>	ELISA, WB, IF, IP, IHC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	HIV-1B

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	10
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein G purification
<b>Dilution Range</b>	ELISA 0.25-0.5 µg/ml; WB 2 µg/ml. Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
<b>Formulation</b>	PBS pH 7.4, with 0.1% sodium azide
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Recombinant Rev protein of HIV-1 BRV (purified from E. coli). Mapped to amino acids 34-40 (TRGARRN – subtype B a.a. consensus)
<b>Immunogen Region</b>	34-40aa
<b>Specificity</b>	Reacts with HIV-1 Rev of subtype B. Mapped to amino acids 34-40 (TRGARRN – subtype B a.a. consensus)
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081