

Human IL2RA (IL2RA) Sandwich ELISA Kit (STJE0000012)

STJE0000012

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL2RA Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	5.89 pg/mL
Detection Limit	31.2-2000pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	IL2RA
Symbol	
Uniprot ID	IL2RA_HUMAN