

Mouse TNF (TNF) Sandwich ELISA Kit (STJE0000017)

STJE0000017

GENERAL INFORMATION

Product Type	Kits
Short Description	This TNF Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	9.56 pg/mL
Detection Limit	31.2-2000 pg/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	TNF
Symbol	
Uniprot ID	TNFA_HUMAN