

Human CSF2 (CSF2) Sandwich ELISA Kit (STJE0000027)

STJE0000027

GENERAL INFORMATION

Product Type	Kits
Description	This CSF2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	6.29 pg/mL
Detection Limit	1000-15.6pg/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	CSF2
Symbol	
Uniprot ID	CSF2_HUMAN