

## Human Tgfb1 (Tgfb1) Sandwich ELISA Kit (STJE0000032)

STJE0000032

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This Tgfb1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| <b>Applications</b>      | ELISA   |
| <b>Reactivity</b>        | Human   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Detection</b>           | Colorimetric  |
| <b>Sensitivity</b>         | 11.21 pg/mL   |
| <b>Detection Limit</b>     | 31.2-2000 pg/ml   |
| <b>Assay Time</b>          | 4.0 hours   |
| <b>Storage Instruction</b> | Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert. |
| <b>Note</b>                | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |             |
|---------------------|-------------|
| <b>Sample Types</b> |             |
| <b>Gene</b>         | TGFB1       |
| <b>Symbol</b>       |             |
| <b>Uniprot ID</b>   | TGFB1_HUMAN |