

Human IL9 (IL9) Sandwich ELISA Kit (STJE0000056)

STJE0000056

GENERAL INFORMATION

Product Type	Kits
Short Description	This IL9 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	8.46pg/mL
Detection Limit	15.6-1000pg/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	IL9
Symbol	
Uniprot ID	IL9_HUMAN