

Porcine CTX1 (Cross Linked C-telopeptide of Type 1 Collagen) Sandwich ELISA Kit (STJE0006677)

STJE0006677

GENERAL INFORMATION

Product Type	Kits
Short Description	This CTX1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in porcine cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Porcine

PRODUCT PROPERTIES

Detection	Colometric
Sensitivity	0.188ng/mL
Detection Limit	0.313~20ng/mL
Assay Time	3.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, plasma and other biological fluids
Gene Symbol	
Uniprot ID	