

Human IgA (Immunoglobulin A) Sandwich ELISA Kit (STJE0006750)

STJE0006750

GENERAL INFORMATION

Product Type	Kits
Short Description	This IgA Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colormetric
Sensitivity	0.14ng/mL
Detection Limit	0.78-50ng/mL
Assay Time	1.5h
Storage	If unopened the kit may be stored at 2-8 °C for up to 6 months. Following first use any unused wells and reagents should be stored according to the components table in the manual.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	Serum, Plasma, Urine, Saliva
Gene	
Symbol	
Uniprot ID	