

Mouse MAdCAM1 (Mucosal Addressin Cell Adhesion Molecule 1) Sandwich ELISA Kit Ready-To-Use (STJE0016849)

STJE0016849

GENERAL INFORMATION

Product Type	Kits
Short Description	This MAdCAM1 Sandwich ELISA Kit, Ready-To-Use is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in mouse cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Mouse

PRODUCT PROPERTIES

Detection	0.156-10ng/mL
Sensitivity	0.059ng/mL
Detection Limit	0.156-10ng/mL
Assay Time	3 hrs
Storage	The whole kit may be stored at -20°C for up to 16 months from receipt. An unopened kit may be stored in the fridge at 2-8°C for up to 6 months. Once opened store individual kit contents according to components table provided with the kit.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, tissue homogenates, urine or other biological fluids.
Gene	Madcam1
Symbol	
Uniprot ID	MADCA_MOUSE